

PFAS BOTTLE KITS

BOTTLE/JAR KITS INCLUDE:

- Labeled bottles and/or jars
- Sampling guidelines
- Chains of Custody (COC)
- Return instructions

Kits are shipped in a cooler suitable for returning samples to our lab for analysis.

KEYS TO SUCCESS RETURNING SAMPLES:

1. Fill out COC completely – provide contact phone/email and PO#
2. Follow sampling guidelines
3. Return samples on ice (no blue ice)
4. Maintain a 4°C during shipping

BOTTLE & JAR SIZES

QTY

• 8 mL*	_____
• 30 mL*	_____
• 60 mL*	_____
• 250 mL	_____
• 4 oz jar	_____
• Trip Blank	_____
• Field Blank	_____
• Extra COCs	_____
• 1L PFAS Free Water (\$25 ea)	_____

* suitable for high level PFAS samples and PFAS screen analysis

SHIPPING: Enthalpy Analytical ground ships bottle kits domestically to our clients free of charge with one business day notice. Same day expedited shipping requests accepted prior to 2 pm are charged to client.

Company Name: _____

Attn: _____

Address: _____

Phone: _____

Email Address: _____

FEDEX # (expedite requests) _____

allows us to provide FEDEX tracking information to you

RETURN SAMPLES TO:

Enthalpy Analytical, LLC, ATTN: Receiving
2714 Exchange Dr., Wilmington, NC 28405



Return completed form to: PFASBottleKit@enthalpy.com